


When do you **need** to add a breast pad?

Breast pads are important when there is pitting oedema (when finger pressure leaves a lingering indentation in the skin) and/or fibrosis (stiffness or scarring) resulting from surgery or radiation. Both conditions impede the natural flow of fluid and contribute to chronic swelling in the breast or chest wall.



Axilla Pad

The **Axilla Pad** was created in response to a request by a breast surgeon looking for a "tapeless solution" immediately following axillary node dissection. The flat pad is placed against the chest wall and held in place by the **Bellisse®**.



Colour



Buff

Fabric

Organic Cotton
& Spandex with
elastic band

The size is determined by the widest biceps measurement:
30-38 cm or 38-46 cm

REF-NO.	Size
75925-02	Axila Pad 30-38 cm Small/ Medium Buff
75925-03	Axila Pad 38-46 cm Large/ Xlarge Buff