Aspire Overbed Tables are designed to meet the needs of users and carers in the Home, Aged Care and Hospital environments.

The robust steel frame and the reinforced base mount provide a very strong, stable base and table top.

For ease of use a quick adjustment lever (spring assisted) and two locking castors are included as standard.

The C-shape base is suitable for use around chairs and under beds. Comes with a beige frame and choice of either Laminate top in White or Beech, or Thermoformed Recessed top in Cream.

## LAMINATE TABLE TOP

accommodates food trays, personal items and provides a stable writing surface

### **FLAT PACK CARTON DESIGN**

is utilised to improve shipping and storage efficiency and reduce potential for damage



Laminate - Beech



Recessed - Cream



RECESSED TOP collects and prevents liquids from spilling over the edge





QUICK
ADJUSTMENT
LEVER makes
height adjustment
easy for users
and carers



**FIXING KNOB** to tighten upright in place



### 12 Month Warranty

Code	Model	Table Top Length	Table Top Width	Height Range
BEA006700	Overbed Table Laminate - AC White	800 mm	430 mm	730 – 1095 mm
BEA006701	Overbed Table Laminate - Beech	800 mm	430 mm	730 – 1095 mm
BEA028000	Overbed Table Recessed - Cream	800 mm	430 mm	730 – 1095 mm



For your nearest Aspire distributor call

1300 133 120

www.aspirecare.com.au



# aspire overbed table SPLIT TILT TOP

The Aspire Overbed Table with Split Tilt Top has a top surface which can be tilted and fixed in place to suit the user while the smaller top surface remains horizontally fixed to hold drinks or other items. This can be useful for writing, drawing or using a tablet device.

The same robust frame and reinforced base mount provide a strong and stable platform while the locking adjustment levers eliminate top movement when locked in place.

### **THERMOFORM SPLIT TILT TOP**

accommodates an easy to clean surface which can be tilted for reading or other activities

### **FLAT PACK CARTON DESIGN**

is utilised to improve shipping and storage efficiency and reduce potential for damage



QUICK **ADJUSTMENT LEVERS** makes height and tilt adjustment easy for users and carers



**FIXING KNOB** to tighten upright in place



Recessed - Cream

**RECESSED TOP** collects and prevents liquids from spilling over the edge

> TWO LOCKING CASTORS allow the table to move easily and be positioned securely where it's needed most



Code	Model	Table Top Length	Table Top Width	Tiit Section Length	Fixed Section Length	Height Range
BEA028100	Overbed Table Split Tilt Top - Recessed - Cream	800 mm	430 mm	600 mm	200 mm	780 - 1140 mm

👣 aspire

For your nearest Aspire distributor call

1300 133 120

www.aspirecare.com.au